

**G<sub>0</sub>/G<sub>1</sub> analysis of KLS cells**  
Bone marrow

1. Remove Bones, clean and crush in ice cold HBSS 2% FBS in a mortar and pestle.
2. Filter cells through 40µm filter to remove bone fragments, pellet cells at 380 x g at 4° C for 4 minutes.
3. Resuspend in 1mL ACK red cell lysis, 5 minutes room temperature, stop with 5mL HBSS+, pellet cells at 380 x g at 4° C for 4 minutes.
4. Resuspend in 5mL HBSS+, count 1:20 in trypan blue.
5. Stain cells using 1µl Lineage mix/ 10<sup>6</sup> cells.

**Lineage Mix**  
(Caltag except IL-7Rα eBiosciences)

CD3	CD4	CD8	GR-1	B220	CD19	IL-7Rα	TER119	MAC-1
1	2	1	2	2	1	2	1	2

6. Wash cells 2x with 1mL HBSS+, pellet cells at 380 x g at 4° C for 4 minutes.
7. Resuspend cells in 2mL HBSS+ in a FACS tube, then add 100µl Dynal beads. Dynal beads should be adjusted according to manufacturers recommendations. Nutate in a covered FACS tube at 4° C for 25 minutes.
8. Deplete lineage<sup>+</sup> cells by incubating FACS tube in Promega MagneSphere® Technology Magnetic Separation Stand (Z5333) for 2 minutes. Remove lineage<sup>neg/low</sup> cells with pasteur pipet.
9. Wash cells with 1mL HBSS+, pellet cells at 380 x g at 4° C for 4 minutes, resuspend in 1 mL HBSS+, count 1:10 in trypan blue.

Expect 10-20% of starting cell number post depletion.

10. Restain cells using 1µl Lineage mix/ 10<sup>6</sup> cells.
11. Wash cells 2x with 1mL HBSS+, pellet cells at 380 x g at 4° C for 4 minutes

12. Stain with  $2\mu\text{l} / 10^7$  cells Tricolor anti Rat (Caltag) for 20 minutes on ice.
13. Wash cells 2x with 1mL HBSS+, pellet cells at  $380 \times g$  at  $4^\circ\text{C}$  for 4 minutes.
14. Resuspend in  $150\mu\text{l}$  HBSS+, add  $50\mu\text{l}$  Rat Ig to block for 5 minutes in the dark.
15. Add  $10\mu\text{l}$  Sca-1/FITC, cKit/APC (BD Pharmingen), CD48/PE (eBioscience) mix /  $10^7$  cells, incubate for 20 minutes on ice in the dark.
16. Wash cells 2x with 1mL HBSS+, pellet cells at  $380 \times g$   $4^\circ\text{C}$  for 4 minutes.
17. Resuspend in 1mL HBSS+ with 1mg/mL

**Sort controls**

$2 \times 10^6$  each sample

1. Unstained
  2. CD44/FITC
  3. B220/APC
  4. FcR/PE
  5. Lin mix+Tricolor
  6. cKit/APC+Sca-1/FITC+CD48/PE+PI
  7. Lin mix+Tricolor+Sca-1/FITC+CD48/PE+PI
  8. Lin mix+Tricolor+cKit/APC+CD48/PE+PI
  9. Lin mix+Tricolor+Sca-1/FITC+cKit/APC+PI
  10. Lin mix+Tricolor+cKit/APC+Sca-1/FITC+CD48/PE+PI
18. Stain sort controls along with sort samples.
  19. Use  $2.5\mu\text{l}$  lineage mix and  $1\mu\text{l}$  other antibodies except  $0.25\mu\text{l}$  CD48 per control.
  20. Sort KSL CD48<sup>+</sup> and CD48<sup>-</sup> cells for each sample.
  21. Spin at 1480 RPM for 5 min ( $441 \times g$ ).
  22. Incubate in  $300\mu\text{l}$  Hoechst  $10\mu\text{g/ml}$  and  $1\mu\text{l/ml}$  (verapamil  $100\text{mM}$  Molecular Probes) in Hoeschst staining buffer for 45 minutes at  $37^\circ\text{C}$ .

23. Add 6 $\mu$ l (Pyronin Y 100 $\mu$ g/ml) incubate for an additional 15 minutes at 37°C.
24. Wash cells 2x with 1mL HBSS+, pellet cells at 441 x g at 4° C for 4 minutes.
25. Resuspend in 300 $\mu$ l HBSS+

**G<sub>0</sub>/G<sub>1</sub> controls**

BM Hoechst + Pyronin Y  
BM Pyronin Y  
BM Hoechst  
BM unstained  
Thymus unstained  
Thymus Pyronin Y  
Thymus Hoechst  
Thymus Hoechst + Pyronin Y

Hoechst staining buffer

45 ml HBSS  
5 ml FBS  
50mg dextrose/ glucose  
1 $\mu$ l 1M HEPES