

Ca²⁺ Measurements in Primary Mouse Lymphocytes

Procedure:

Design experiment

- Decide which tissues and cell types to analyze in how many mice
 - ==> Calculate sample numbers and prepare sufficient reagents and solutions
- Include the following compensation controls to set compensations and optimize PMV voltage settings:
 - unstained cells
 - PE AB only
 - APC AB only
 - {Third color} AB only
 - Indo-1 (Ca²⁺ dye mix) only

Collect baseline data in MMF first, then add Ca²⁺-add back and Ionomycin to collect data for maximal Ca²⁺ influx. Adjust PMV voltage to optimally separate blue and violet Indo-signal. Then compensate the other channels under maximal Ca²⁺ influx.

Prepare lymphocyte suspensions

- conduct this part at RT
- harvest desired organs (spleen, LN, BM, thymus) in 5 ml M199-2 (10mM HEPES, 2% FCS, 1X Pen/Strep)
- gently disperse cells in M199-2 through a cell strainer with the plunger of a 1 or 3 ml syringe into a 50 ml tube
- wash strainer with 1ml M199-2 and add the remaining cells to the tube
- spin 5 min at 1500 RPM, discard supernatant
 - a) If working with spleen or BM, or if significant contamination with red blood cells occurs:
 - resuspend cells gently in 2 ml 1x BD Pharmlyse (cat # 555899); incubate at RT for 2-3 min.
 - add 10 ml M199-2 and mix gently
 - immediately spin 5 min at 1500 RPM, discard supernatant
 - resuspend cells in 5 ml M199-2
 - b) If no RBC contamination occurs, then proceed to cell counting
- count cells
- spin 5 min at 1500 RPM, discard supernatant
- resuspend cells required for assay in MFC at 10⁷ cells/ml; for remaining cells, keep in M199-2
- ==> At this point, cells can be stored at RT for a few hours (but not > 6 hrs)

Prepare Ca²⁺ dye solutions

- for 2µg/µl (=2 mM) **Indo-1 stock**, reconstitute 50 µg Indo-1 in 25 µl DMSO at RT, vortex, spin briefly.
- you will need 5 µl stock solution per sample
- keep at RT in the dark (drawer, closet)

- make fresh **2x CDLS** from this stock for each sample batch (2.5 µl Indo-1 stock in 1 ml MC, use within 1 hr). You will need 250 µl of this solution per sample.

Loading the cells with Ca²⁺ dye and staining for surface markers

Ca²⁺ flux profiles may vary significantly between cells that are run immediately after loading and cells from the same batch that have been loaded >2 hrs before. Therefore, never load more samples than you can run within 2hrs (typically 4-6 samples plus loading controls)

- For the following, use 2.5 x 10⁶ cells (250 µl of the above cell stock) per sample. Samples that will be stained for the same surface markers can be pooled.

- make the 2x CDLS fresh for each sample batch (One freshly made dye stock solution in DMSO can be used for all samples on the same day)

- transfer 2.5 x 10⁶ cells into individual microtubes

- add 1 vol of CDLS

- incubate for 30-45 min at 37 °C in the dark (bacterial incubator). Invert gently to mix every ~10 min.

- *During this incubation time, prepare surface marker staining and stimulation mixes (see below)*

- spin loaded cells for 5 min at 1500 RPM at RT, discard supernatant

- add 900 µl MFC to loaded cells to wash

- repeat wash (MFC) and resuspend in 200 µl MFC

- *stain the cells on ice for surface markers, use cold buffers and treat very gently to preserve membrane integrity*

- add 1/10 vol (25 µl) FACS-AB premix

- **Stimulation variant A:** Add 1/10 vol (25 µl) Stimulation-AB Mix1A here

(In the alternativestimulation variant B, this is instead done during FACS analysis)

- incubate 15 min on ice in the dark

- *during this time, go to Flow facility, turn on 37 °C waterbath and UV laser on FACS machine*

- wash cells 1-2x in MFC

- wash cells 1x in MMF

- resuspend in 500 µl MMF and transfer into FACS tube, store on ice in the dark

==> the cells are now ready for FACS analysis in a Ca²⁺ re-addition protocol

FACS analysis

(~15 min per sample, stimulation mixes must be pre-warmed to 37 °C)

- Take the following to the FACS facility (unless already available there):

- Heat block (37 °C)

- Stimulation mixes

Variant A: Mock, SA, Tg, Ca²⁺ add-back, Ionomycin

Variant B: Mock, Stimulation-AB Mix 1B, xABMix, Tg, Ca²⁺ add-back, Ionomycin


- Cells on ice in the dark (Aluminum foil on top of ice bucket)

- Pipets (P200, P1000), yellow & blue tips or better filter tips

- 10% bleach
- deionized H₂O

- Check/adjust filter settings (Here shown for LSRII):

- FSC: Linear
- SSC: Linear
- PE: Logarithmic
- APC: Logarithmic
- Indo-1 violet: Linear
- Indo-1 blue: Linear



- Set compensations (not needed for combination PE/APC/Indo)
- Read Violet/Blue ratio over time (and plot for monitoring progress)

- place one cell sample at 37 °C for 5 min (during which you bleach, prime and flush with water)
- run bleach for 2 min, fast speed
- run water for 2 min, fast speed
- prime
- set speed to low
- make sure FACS software is set to new sample with correct annotation
- collect **30 sec** (*baseline*)

- **Stimulation variant B:** *This is not done if stimulation variant A is used!*

- add 100 µl Stimulation-AB Mix 1B (5 µg/ml final in MMF), vortex or swirl briefly
- collect **1 min** (*AB baseline*).
- add 50 µl Stimulation Mix 2 (Mock or SA or xAB or Tg), vortex or swirl briefly
(*If machine runs out of sample before time course over, change SA addition to adding 400 µl of 1:10 diluted Stimulation Mix 2, which gives ~ the same final conc in a large enough volume for a 15 min run*)
- collect **3 min** (*intracellular store release*)
- add 50 µl Ca²⁺ add-back, vortex or swirl briefly
- collect **5 min** (*plasma membrane influx*)
- **If SOC/CRAC modulating compounds are being tested**
 - collect only **3 min** (*plasma membrane influx*)
 - add 50 µl mock or compound at desired concentration in MFC, vortex or swirl briefly
 - collect **2 min** (*modulated plasma membrane influx*)
- add 50 µl Ionomycin
- collect **1-2 min** until peak flux reached (*loading & viability control*)
- place next cell sample at 37 °C for 5 min and start again by bleaching, flushing and

- If last sample done: Bleach, flush, prime and then place sample capillary into water, set machine on standby.

- Export FCS files onto server, memory stick or CD/DVD and analyze offline in Flowjo.

- If more samples will be run on same day, keep machine on standby. Otherwise, turn off machine.

Media and solutions: (FCS can be replaced with FBS)

CM: RPMI + 10 mM HEPES + 10% FCS + 1% Glu (+ Pen/Strep/Glu = PSG)

M199: M199 + EBSS + 20 mM HEPES + 2% FCS + 1% Glu + 2.2 g/L Na-Bicarbonate (+ PSG)

MFC: HBSS with Ca^{2+} and Mg^{2+} (~ 1 mM each) + 10 mM HEPES + 1% FCS (phenol-red free)

MC: MFC without FCS

MMF: HBSS without Ca^{2+} and Mg^{2+} + 10 mM HEPES + 1% FCS + 1 mM MgCl_2 + 1 mM EGTA (phenol-red free)

Indo-1 stock, 2 $\mu\text{g}/\mu\text{l}$ (=2 mM): Reconstitute 50 μg Indo-1 in 25 μl DMSO at RT, vortex, spin briefly, store at RT in the dark.

2x Ca^{2+} dye loading solution (2x CDLS): 2.5 μl Indo-1 stock in 1 ml MC, use within 1 hr, store at RT in the dark. The final loading conc will be 5 $\mu\text{g}/\text{ml}$. *If loading is not good, try 10 $\mu\text{g}/\text{ml}$.*

10x FACS AB mixes:(need 25 μl per sample)

- T cell mix: $\alpha\text{CD4-APC}$, clone GK1.5, [1:25] + $\alpha\text{CD8-PE}$, clone 53-6.7, [1:25] in MFC
- B cell mix: $\alpha\text{CD4-APC}$, clone GK1.5, [1:25] + $\alpha\text{CD8-APC}$, clone 53-6.7, [1:25] + $\alpha\text{B220-PE}$ (or CD90-PE) [try 1:25] in MFC
- B/T cell 5-color stain: $\alpha\text{CD4-APC}$, clone GK1.5, [1:25] + $\alpha\text{CD8-PE}$, clone 53-6.7, [1:25] + αB220 or αCD90 -*{determine suitable color, must be compatible with PE, APC, Indo-blue and Indo-violet}* [try 1:25] in MFC. This can then replace (a) and (b).

10x Stimulation-AB Mix 1A:(need 25 μl per sample)

- CD3 stimulation, mild: $\alpha\text{CD3-Bio}$, clone 2C11, [1:36] in MFC
- CD3 stimulation, strong: $\alpha\text{CD3-Bio}$, clone 2C11, [1:4] in MFC
- CD3/CD4 co-stimulation: $\alpha\text{CD3-Bio}$, clone 2C11, [1:4] + $\alpha\text{CD4-Bio}$, clone RM4-4, [1:10] in MFC
- BCR stimulation: $\alpha\text{IgM-Bio}$, [1:xxx] or $\alpha\text{IgD-Bio}$, [1:xxx] in MFC *{Optimize, start at 1:4}*
F(ab')₂ α -mouse IgM at a final concentration of 30 $\mu\text{g}/\text{ml}$ in MFC.

10x Stimulation-AB Mix 1B:(need 100 μl per sample)

- CD3 stimulation, mild: $\alpha\text{CD3-Bio}$, clone 2C11, [1:xxx] in MMF
{Optimize - try 1:36 to start}
- CD3 stimulation, strong: $\alpha\text{CD3-Bio}$, clone 2C11, [1: xxx] in MMF
{Optimize - try 1:4 to start }
- CD3/CD4 co-stimulation: $\alpha\text{CD3-Bio}$, clone 2C11, [1: xxx] + $\alpha\text{CD4-Bio}$, clone RM4-4, [1: xxx] in MMF
{Optimize - try 1:4 respectively 1:10 to start}
- If unbiotinylated, purified ABs are used: Make a stock in MMF that yields 5 $\mu\text{g}/\text{ml}$ final
{Optimize}
- BCR stimulation: $\alpha\text{IgM-Bio}$, [1:xxx] or $\alpha\text{IgD-Bio}$, [1:xxx] in MMF*{Optimize, start at 1:4}*

10x Stimulation Mix 2:

Mock: MMF

10x SA to crosslink Bio-ABs: 50 $\mu\text{g}/\text{ml}$ streptavidin in MMF, need 50 $\mu\text{l}/\text{sample}$, final conc in tube ~5 $\mu\text{g}/\text{ml}$). SA stock is 10 mg/ml.

xAB mix to crosslink purified non-biotinylated ABs: Make an anti-isotype AB stock in MMF that yields 5 $\mu\text{g}/\text{ml}$ final

10x Tg: 20 μ M Thapsigargin in MMF, need 50 μ l/sample (final conc in tube \sim 2 μ M)

10x Ca²⁺ add-back: 20 mM Ca²⁺ in MFC (will reconstitute \sim 1 mM Ca²⁺ due to presence of EGTA in MMF. 20 mM Ca²⁺ can be made by diluting 40 μ l 1 M CaCl₂ into 2 ml MFC. If it doesn't work well, try 50 mM Ca²⁺ in MFC.)

Ca²⁺ Influx modulating compound: Add xy x desired final concentration in MFC. *Determine factor xy by measuring remaining cell volume in a mock Flow experiment!*

Tentatively, try 10x cpd in MFC (for example 5 mM or 1 mM IP₄/AM or IP₄/PM)

lonomycin: 50 μ g/ml Ionomycin in MFC, need 50 μ l/sample (final conc in tube \sim 5 μ g/ml)

Reagents

Bio-CD3 BD#553060 (2C11), 0.5 mg/ml
Bio-CD4 BD#557443 (RM4-4), 0.5 mg/ml
CD4-APC Other epitope for FACS: eBiosc #17-0041-83 (GK1.5) 0.2 mg/ml
CD8a-PE eBiosc #12-0081-83 (53-6.7) 0.2 mg/ml
B220-xxx xxx
F(ab')₂ α -mouse IgM Jackson Immuno Research xxx (?)

α IgM-Bio xxx
 α IgD-Bio xxx

Ca²⁺ dyes: (Make in MC buffer)

Indo-1 Molecular Probes #11223, 20 x 50 μ g

OR

Fura Red, 2 mL/mL (1mM stock, Molecular Probes Cat: F-3021)

Fluo 4, 1 mL/mL (1mM stock, Molecular Probes Cat: F-14201)

Pluronic, 1 mL/mL (20% Solution, Molecular Probes Cat: F-127)

Use Fluo-4 at 3 mg/ml and Fura Red at 5 mg/ml, load in RPMI/3% FCS

SOC Modulators:

Tg: Thapsigargin, xxxx, Sigma (?)

Ionomycin: Ionomycin, xxxx, Sigma (?)

Exemplary reagent addition scheme for a Ca²⁺ add-back experiment

(Here, 400 μ l Stimulation Mix 2 are added and the times are slightly different from above. They can be adjusted as needed.)

