

Cystogram

Scout radiograph	Images	Purpose	Example
Yes	AP supine pelvis	Foley position, post op changes, extraluminal air, fractures (if trauma)	
Positions and contrast	Images	Purpose	Example
Water soluble contrast Fill bladder via gravity drip until patient discomfort or flow stops	Fluoro store as bladder fills	Look for early extravasation Look for ureteroceles or other filling defects	
	AP Radiograph (using the II)	Look for extravasation	

	LPO Radiograph (using II)	Look for extravasation	
	RPO Radiograph (using II)	Look for extravasation	

Drain contrast	AP Post Void Radiograph (overhead or II)	? extravasation that was hidden behind full bladder	
Optional	Images	Purpose	Example
	AP Fluoro store or Radiograph (II) kidneys	For vesicoureteral reflux	